




Alger Equipment Co., Inc.  
320 Flightline  
Lago Vista, TX 78645


Ph 512.267.0383  
800.320.1043  
Fax 512.267.0452

**Alger Light**  
with **LED**



### Alger Light Positioning (Side View)

Ideally position the Alger Light with LED between the hip and knee if you work with the patient in the "reclined" position.



### Alger Light Positioning (Overhead View)

Ideally center the Alger Light with LED over the center-line of chair. If there is obstruction in the center the light can be off-center approximately 12" to either side of center.

